CALIBRADOR DE COAGULACIÓN

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COAGULATION CALIBRATOR

COD 61006 4 x 1 mL
STORE AT 2-8°C
Only for <i>in vitro</i> use in the clinical laboratory

PRODUCT DESCRIPTION

The Coagulation Calibrator is a lyophilized pooled human plasma containing component concentrations suitable for the calibration of measurement procedures.

CONTENTS AND COMPOSITION

Coagulation Calibrator. For 1.0 mL. Lyophilized human plasma.

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

INSTRUCTIONS FOR USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 1.00 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 30 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. The reconstituted material is to be treated like the patient plasma.

STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Coagulation Calibrator is stable until the expiration date given in the label.

The components of the reconstituted Coagulation Calibrator are stable for 4 hours at 22°C.

Reconstituted Coagulation Calibrator is stable for at least 8 days at -20°C. Freeze only once.

SPECIFIED VALUES

The assigned concentration values for components and their traceability are shown in the enclosed tables. Traceability can be assured only if the BioSystems reagents and recommended measurement procedures are used.

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